

BEAMER

For advanced electron beam lithography

Marvin Zai, PhD
zai@genisys-gmbh.com

- GenISys Introduction
- E-Beam Lithography: From Layout to Wafer
- BEAMER Basics
- Basic Pattern Preparation
- Fracturing
- Field Control
- Proximity Effect Correction (PEC)

GenISys offers software solutions for optimization of micro- and nano-fabrication processes

Company:

- Founded in 2005, joined RSBG Group in 2018
- Headquartered in Taufkirchen - Munich, Germany
 - Additional development location in Jena – Germany
 - Subsidiaries for customer support in USA and Japan
- Fast, Flexible, Responsive



Foto Sterflinger

Electron and Laser Beam Direct Write Software

- Market leader for Gaussian beam direct write systems
- Installed at most major nano-fabrication centers worldwide, has become a MUST for advanced e-beam lithography



Monte Carlo simulation software

- MC-Simulation of electron distribution for e-beam lithography modeling and correction
- Process Calibration, PSF visualization, extraction and management



3D lithography simulation & OPC software

- Proximity Lithography (mask aligner) & Projection Lithography (stepper / scanner)
- Electron Beam Lithography, Laser Beam Lithography (Heidelberg Instruments laser systems)



SEM Image Analysis & Metrology

- Metrology software for SEM based metrology and inspection



Mask Production Software

- Dedicated MDP for mask house, high performance (hierarchy, parallel processing, mask process correction...)
- Special Application: Flat Panel Display, Photonic IC, non-IC



Selected Installed Base

- Major nanofabrication centers worldwide
 - Universities, Research Centers
- Industrial R&D and special production
 - Advanced FPD manufacturers
 - Mask manufacturer



500 commercial licenses worldwide

- 230 BEAMER (83 EU, 80 APAC, 70 US)
- 140 TRACER
- 65 LAB
- 30 ProSEM

Strategic Partnerships

GenISys is an independent software supplier working with all major lithography and inspection system manufacturers.

RAITH
NANO FABRICATION

JEOL

ELIONIX

ADVANTEST

CRESTED

HEIDELBERG
INSTRUMENTS
The power of direct writing

SÜSS + MicroOptics


E V G


ZEISS





KBTEM-OMO

- GenISys Introduction
- E-Beam Lithography: From Layout to Wafer
- BEAMER Basics
- Basic Pattern Preparation
- Fracturing
- Field Control
- Proximity Effect Correction (PEC)

Electron Beam Lithography

E-Beam lithography (EBL) is the most utilized technology for patterning nano-scale (Quantum) devices

- Beam size down to few nano-meter
- Most flexible pattern and substrate
- Direct Write from CAD data to sample



Source: AMO GmbH - Germany


Source: NIST CNST - USA

Source: Weizmann Institute – Israel
Stanford University, USA

GOOD DATA IN




GOOD SAMPLE OUT




- Electron beam systems are not “office printers”
 - Many application and process specific effects to fight
- Easiest (but expensive) way until now: trial-and-error
 - Using smallest possible beam current, accepting long write times
 - Required many iterations to determine best exposure dose
- BEAMER established new data preparation “best practices”
 - Enables faster and more accurate writing using advanced error correction and data preparation techniques
- Good data preparation requires understanding key concepts...

- The data (GDSII, DXF, etc.) is your design, the layout.
- The act of converting the data to the machine format is often referred to as fracturing.
- What is fracturing?
 - Fracturing is the method by which a complex shape is broken down into simple (primitive) shapes (trapezoids).
 - Most e-beam tools can only accept trapezoids





Shape Fracture/Exposure



Exposure in Write Fields


1.62 x 1600 nm
Exposure requires scanning




the write sequence of field by field

mechanical)
on (electro-static or –
n has limited range
en accuracy and
e move / no scan is slow
/ stage steps is faster
es you control over

Electron Scattering




E-beam Lithography Distortion



What you get is
NOT what you
design

Result on sample after etch





BEAMER enables clever data preparation...

- Mitigate e-beam artifacts
 - Shot Filling, Field Positioning Errors (aberration/stitching)
- Correct for proximity effects
 - Long-Range, Short Range, Process Effects, 3D PEC
- Improve
 - Wafer: Resolution, Linearity, Line Edge Roughness, Uniformity
 - Process: Repeatability, Write Time
- Essential for writing high precision devices

- GenISys Introduction
- E-Beam Lithography: From Layout to Wafer
- BEAMER Basics
- Basic Pattern Preparation
- Fracturing
- Field Control
- Proximity Effect Correction (PEC)


Unique VisualFLOW user interface

- Comprehensive functional library
- Easy and fast operation
- Supports Windows & Linux
- Flexible licensing



Supports all major data formats


- Import / Export of all major layout formats
 - GDSII
 - OASIS
 - DXF (AutoCad)
 - CIF and Ltxt (text-based formats)
 - Bitmap
- Export / Import of all major e-beam machine formats (optional)
 - Raith EBPG (GPF/IWFL), VB (VEP/FRE), Voyager, Raith GDS
 - JEOL (v30, j51, j01)
 - Elionix, Crestec (cel, con)
 - Heidelberg Laserwriter
 - ADVANTEST e-beam lithography system
 - MEBES mask writer
 - KBTEM exposure systems



BEAMER offers a rich library of data preparation functionality

- Extract layout elements, layers, cells, regions
- Filter layout elements by width, height, area
- Transform layouts: scale, mirror, rotate, shift
- Grid optimization control
- Map layer names / numbers
- Fracture with extensive controls (machine independent)
- Heal layout elements or remove overlaps
- Reverse tone
- Bias or resize intelligently
- Apply Boolean operations: XOR, P-XOR, OR, AND, MINUS
- Merge layouts


Correction and Simulation




BEAMER offers comprehensive correction and simulation:

- Standard Dose Proximity Effect Correction (PEC)
- Shape PEC
- 3D PEC
- Corner PEC
- FDA
- Rule-OPC
- E-Beam simulation
- Metrology


- We call each BEAMER “project” a flow
- Flows are constructed of connected modules
- A module takes an input pattern and produces a single output pattern
- The data “flows” from module to module along each connector line
- The processed pattern is exported to a chosen file format




The status of a module is shown by the module color and the icon marker located at the right side of the module.

Module State	Description
 Import	A module which has not had its parameters set is identified as a module missing its inputs and outputs.
 acmos	The minus sign to the right of the module signifies that parameters have been set and a layout named acmos is chosen but the module has not been run yet.
 acmos	A progress bar is shown while the process is being computed and the time passed is shown on the bar to the right.
 acmos	A module which has been run and the results are ready, is shown with a box View icon.


Graphical User Interface



- Sample flows, sample patterns and training flows in Examples tab
- Drag and drop into flow
- Every flow is pre-set with all needed parameters and patterns




- Help Manual, Release Notes and Formatter Specific Notes available under “Help”
- Select the needed formats to send PDF versions of these to a directory for later viewing.




Customizable Defaults

- Choose the settings item on the right side of BEAMER to reach customizing properties.
- Resetting defaults can also be done through the context-menu.





- GenISys Introduction
- E-Beam Lithography: From Layout to Wafer
- BEAMER Basics
- Basic Pattern Preparation
- Fracturing
- Field Control
- Proximity Effect Correction (PEC)




- Heal removes overlap between shapes
- Two processing modes:
 - “Healing” removes overlap and merges adjacent shapes
 - “Overlap Removal” does not merge adjacent shapes
- Can be applied on all layers, individual layers or by dose





- Correct for process bias by “growing” or “shrinks” shapes by adjusting shape edges
- All layers or individual layers
- Be careful...
 - A negative bias can erase parts of a layout.
 - Shapes that “grow” into each other are merged
 - Adjacent shapes are not separated




- Resizes shapes by moving every edge like bias
- Unlike BIAS:
 - Does not merge shapes that “grow” into each other
 - Will pull apart adjacent shapes




NOT allows tone reversal of one input pattern.




AND keeps what's in common between two inputs patterns.



OR performs a merge and heal between shapes with two input patterns.






XOR removes any areas that are over-layed on the two input patterns.



A	B	Q
0	0	0
0	1	1
1	0	1
1	1	0


- MINUS subtracts input B (on the right) from input A (on the left)
- Order matters.






- Odd overlays (1,3,5, ...) are kept.
- Even overlays (2,4,6, ...) are removed.

Correction and Simulation




Transform allows standard pattern transformations

- Center
- Rotate
- Reflect
- Scale
- Shift



Extract allows shape extraction while preserving hierarchy




- Cell Extraction
 - Instances
 - Definitions
- Layer Extraction
- Region Extraction
 - By Region Layer
 - Visually
 - Start/End Region Selection: Shift + Left Mouse Click
 - Region coordinates are placed in the table to the left




Region Extract Modes

Region Extract Options:


- Clip cuts all elements at the region border
- Within extracts only polygons that are completely inside the region
- Touching extracts all polygons that are partially inside the region



Clip




Within




Touching

- Changing “Database Grid” adjust the layout’s base unit
 - Vertices must be on grid
 - Setting a coarser grid may cause grid snapping
- Merging patterns requires setting the same database grid







1nm grid

No smoothing tolerance



1nm grid

3nm smoothing tolerance



- Layout Smoothing reduces vertex count
 - Smooths shape edges, tolerance is allowed deviation of edges of smoothed shape from designed shape
 - Can reduce fracturing time, post-fracture shape count and total data volume
 - Inspect result!

Opens fully functional Layout Editor


- Edit a pattern in a workflow
- Create a new pattern

Things to remember:

- Saving the Edit module in the flow does not save your pattern if you close and re-open your flow.
- EXPORT and then IMPORT to save modifications or newly created patterns.




Cooperation with
Juspertor



- GenISys Introduction
- E-Beam Lithography: From Layout to Wafer
- BEAMER Basics
- Basic Pattern Preparation
- Fracturing
- Field Control
- Proximity Effect Correction (PEC)



- Fracturing is breaking up single complex polygons into simple “primitive” shapes, usually trapezoids
- E-beam tools expose primitives shape-by-shape
- BEAMER allows controlling:
 - Fractured shape count
 - Optimizing shot filling of fractured shapes
 - Writing order from shape to shape



Fracturing Modes


Three Fracturing
Modes:

- Conventional
- LRFT – Large Rectangle Fine Trapezoid
- Curved



Conventional

- Mostly horizontal or vertical slices of polygons
- Can cause result in some asymmetry and/or shot filling problems




Fracturing Mode: Conventional

Fracturing Mode: LRFT

LRFT


- Fractures with as many Large Rectangles as possible
- Fills rest of shape with Fine Trapezoids
- Can improve symmetry and shot filling





Curved Fracturing

Curved fracturing

- Remove and/or shift shape vertices while considering resolution and beam step size
- Fewer vertices create fewer fractures
- Just enough vertices to create smooth curve is ideal
- Can improve shot filling



Curved Fracturing Tolerance




- Curved Fracturing Tolerance is manually adjusted
- “Allowed deviation” between design and fitted curve
- Inspect result to find optimal value

- GenISys Introduction
- E-Beam Lithography: From Layout to Wafer
- BEAMER Basics
- Basic Pattern Preparation
- Fracturing
- Field Control
- Proximity Effect Correction (PEC)

- A Gaussian beam fills the pattern with discrete pixels. The exposure quality depends on:
 - Proper shape filling
 - Position of shapes in field/sub-field (aberrations/distortions)
 - Accurate field stitching
- Poorly controlling these things can result in:
 - Defective shapes
 - Line edge roughness
 - Shape distortions at field borders
 - Field/sub-field stitching errors


E-Beam Exposure Fields

- Field size is tool limited so large layouts are tiled
- Default field sorting:
 - Uniform tile size
 - Start in one corner and meander fill
- Not always optimal...




Field Placement Types

- Fixed
 - Standard Fixed Grid
- Floating
 - Fields positioned to enclose shapes
- Manual
 - Free positioning of fields




Field Stitching:

- Moving between fields requires a mechanical stage move.
- Elements crossing the field border are split and exposed in two fields.
- Lower quality at field border due to aberrations/distortions.
 - Field center has highest accuracy; field edges have lower accuracy.
- More time/stage movements between feature writing



Issues at Field Borders




Courtesy: Nikolai Klimov, NIST

Issues at Field Borders

Stitching errors more likely if there is more time between writing both parts of a feature

- Tool drift, many stage movements...



Courtesy: Nikolai Klimov, NIST

Field Ordering determines stage motion and impacts delay between fields


Meander X, Meander Y, Raster X and Raster Y

Example:

- X-Meander is not optimal for a grating oriented in Y-direction

Meander X


Raster Y





Field control – Overlap

- The Field Overlap defines “mutual” area between fields:
 - Elements in overlap are not split, but exposed either completely in one or the other
 - Only elements exceeding the overlap area to both sides need to be cut



Standard Fracture
Shapes Cut





Overlapping Fields
No Shapes Cut




Overlapping Fields



Overlapping Fields





- What is best field sorting method for this kind of pattern?
 - Sparse features
 - “Large shapes” (relative to field size)
- Fixed Fields sorting is not optimal...
 - Shapes are not centered in field where write quality is higher
 - Shapes cross field boundaries which can introduce stitching errors




Floating Fields


- Automatic Floating Field sorting is available
- Fields are automatically placed to surround shapes
 - Features are centered in fields if possible
 - Stage moves are optimized





Manual Field Placement

- Manually select fields
 - Region selection as in Extract
 - Use “Region Layer” definitions to select region
- Various Sorting Options
 - Selected Fields Only
 - Selected Then Remainder
 - Remainder With Selected
 - Selected Then Floating










Selected Then
Remainder


Selected fields are written first, then the rest of the pattern is exposed



Remainder with
Selected



The standard field traversal picks up the selected floating fields during exposure.

Manual Field Control: Layers




Multipass works by exposing the same shape multiple times at a fractional dose

- Reduces LER, improves resolution
- Mitigates stitch errors (with offset)



Single pass exposure with stitch errors




Same pattern fractured with multi-pass



Multipass: Theory of Operation

- Writing Errors are either systematic or statistical...
- Statistical errors include:
 - Beam Current Fluctuations
 - Beam Jitter, Beam Drift
 - Stage Position Errors, mechanical vibrations, ...
- Systematic errors include
 - Field Distortion and Field Aberrations
 - Scan Non-Linearities
 - X/Y asymmetries due to discrete spots in Y / dragging in X
 - ...
- Multipass reduces statistical errors by averaging
- In addition, systematic errors can be reduced by offset strategies



- Resist edge at dose to clear threshold
- Accurate dose and shot placement results in correct CD




- Errors can result in...
 - misplaced shots that shift the feature edges
 - fluctuations in dose applied at each shot that causes oversized or undersized features



- Multipass works by exposing the same spot multiple times at a fractional dose
 - (e.g. 4 pass at $\frac{1}{4}$ dose per shot)
- Doses from each pass sum



4X Pass with $\frac{1}{2}$ Field Shift in X/Y






One E-Beam Exposure Field


Left half of
Yellow area
exposed

4th time


- In multi-pass writing, shapes are written multiple times, averaging both temporal and spatial variations, improving uniformity.
- (Use half the dose for each pass)




High LER



Low LER



Medium LER



Single Pass $100\mu\text{C}/\text{cm}^2$

4-pass with stationary
field $100\mu\text{C}/\text{cm}^2$

4-pass with field
shifts $100\mu\text{C}/\text{cm}^2$

Multipass Mode

- Multipass Mode

- 1 – writing with a single pass
- 2 – every element is written twice
- 4 – every element is written quadruple times
- Dose Selective – For each defined dose value the number of passes can be specified

- Mainfield offset



- Defines the shift of a field against another (Units are field sizes)

- Subfield offset

- Defines the shift of a subfield against another (Units are field sizes)

- Layer for Multipass

- Applies multipass only to specified fields



Dose larger	Number Passes
1.0	1
1.5	2
2	3
4	4
6	8

Feature Sorting in Field

- Feature Order is the order shapes are exposed within each field
- No Compaction – Default, takes data as is
- Array Compaction – Compacts data by recognizing identical shapes
- Writing Order – “Nearest neighbor” exposure
- Follow Geometry – Follows long lines, good for contiguous structures
- Random – Random shots, causes sub-field jumping but can reduce charging effects

No Sorting




Follow
Geometry



- GenISys Introduction
- E-Beam Lithography: From Layout to Wafer
- BEAMER Basics
- Basic Pattern Preparation
- Fracturing
- Field Control
- Proximity Effect Correction (PEC)


Electron-Solid Interactions



- Short-Range Effects(5-50 nm)
 - E-Beam: beam blur + forward-scattering
- Mid-Range Effects(30-500 nm)
 - E-Beam: mid-range scatter for complex material
 - Resist blur/diffusion 30-100 nm
 - Development/Etch loading 200-500 nm
- Long-Range Effects(10-35 μ m)
 - E-Beam: back-scattering



Effects can be summarized by a Point Spread Function – PSF:
Spread of Energy from the incidence of exposure.
The adjacent features interact with each other: Proximity Effect

Point Spread Function - PSF



$$PSF(r) = \frac{1}{\pi(1 + \eta)} \left[\frac{1}{\alpha^2} e^{-\frac{r^2}{\alpha^2}} + \frac{\eta}{\beta^2} e^{-\frac{r^2}{\beta^2}} \right]$$


TRACER calculates the absorbed energy distribution over resist thickness and distance.



Calculation of Absorbed Energy

Knowing the PSF, the absorbed energy at any position x can be calculated:


$$E(x) = P(x) \otimes PSF$$



Calculation of Absorbed Energy

Knowing the PSF, the absorbed energy at any position x can be calculated:


$$E(x) = P(x) \otimes PSF$$



Calculation of Absorbed Energy


Knowing the PSF, the absorbed energy at any position x can be calculated:

$$E(x) = P(x) \otimes PSF$$



Absorbed Energy after PEC


PEC corrects so that absorbed energy crosses threshold at pattern edge




Method resolves lines and gaps symmetrically by reducing exposure dose in large areas to reduce the scattering into gaps

Proximity Effect Correction

- Target of correction:
 - Adjust all feature edges to the same absorbed energy: Threshold Dose of resist
- Inside the feature $E(x) >$ Threshold Dose:
 - Resist will remain (negative)
- Outside the feature $E(x) <$ Threshold Dose:
 - Resist will be removed (negative)
- Correction Equation:
 - $E(\text{edge}) = 0.5 = D(x) \otimes \text{PSF}$
- There is large number of PEC algorithms. The strongest algorithms are based on:




Modeling Absorbed Energy




Absorbed Energy in Resist



Assume resist edge falls
at constant threshold
energy




Assume resist edge falls
at constant threshold
energy






Absorbed Energy in Resist




Corrected Exposure Dose









Uncorrected Pattern




Corrected Pattern





- Always give PEC healed data – no overlaps are allowed for proper correction.
- Always apply PEC after any other layout manipulations (e.g. BIAS, TRANSFORM, MERGE...). You want PEC to correct the actual to-be-exposed pattern.
- Choose the appropriate PSF by selecting from the archive, a gaussian approximation or importing a PSF.
- For most cases, that's all you need to do!




- PEC dose class accuracy is selected on the Accuracy tab
- Fractures are only made when it is required to achieve the specified accuracy
- Smaller value = higher accuracy
- Separation between dose factors is approximately accuracy value
- BEAMER will report the achieved dose accuracy



Larger accuracy tolerance reduces shape count




- Isodose Grid sets “fracture grid” that fractures can occur on
- PEC steps from grid point to grid point
 - The value should equal a multiple of the shot pitch



Short Range Correction

- Include Short Range Correction toggles short range PEC
 - Only necessary when features are on the scale of the blur
 - 20-100 nm for Silicon
- Increases:
 - Calculation time
 - Shape count
 - Dose factor “spread”



Effective Short Range Blur

Proximity Effect Correction

General

Accuracy

Advanced

Comment

Correction Layer Selection
Layer(s) *

PSF Representation
Archive Gaussian Approximation Numerical PSF

`g; Substrate: Si; Layers: ; Resist: PMMA 100 nm; Energy [kV]: 100;`

Archive... Global Archive... View Comment...

Effective Short Range Blur FWHM[um]
0.01


Add Gamma [um] 1
 Include Short Range Correction
 Lateral De ...

10nm Blur

Show Energy Density
 Show Cumulative Radial Energy
Behaviour X-Axis: Logarithmic Quadratic Linear
Behaviour Y-Axis: Logarithmic Linear

Cumulative [1 = 100%]


Short Range
Long Range






Separation at 0.1455 um.
No additional Separation necessary.
Additional Information:
Min. layout independent LR dose factor = 0.7314

OK Cancel

Effective Short Range Blur



- FDA manually assigns or changes dose assignments
 - Assign a new dose, multiply existing doses by a factor, or add a value to existing doses.
 - Assign by layer or within defined regions.



Thank You!

support@genisys-gmbh.com



Headquarters

GenISys GmbH
Eschenstr. 66
D-82024 Taufkirchen (Munich)
GERMANY
📞 +49-(0)89-3309197-60
📠 +49-(0)89-3309197-61
✉️ info@genisys-gmbh.com

USA Office

GenISys Inc.
P.O. Box 410956
San Francisco, CA
94141-0956
USA
📞 +1 (408) 353-3951
✉️ usa@genisys-gmbh.com

Japan / Asia Pacific Office

GenISys K.K.
German Industry Park
1-18-2 Hakusan Midori-ku
Yokohama 226-0006
JAPAN
📞 +81 (0)45-530-3306
📠 +81 (0)45-532-6933
✉️ apsales@genisys-gmbh.com